


ANALYTICAL BALANCE

0.001g (1mg)

TOPLAB INDIA
"Serving Science with Solution"

The visual appearance of the product is blent into the simple concept, stylish and generous. Specially makes the fine pouring cast- aluminium outer covering, so as to enhance the balance of the anti-static anti-jamming capability. The panel with five operational button, Sliding glass windscreen provide big space, insightful vision and operate easily. We use high precision magnet sensor to keep the weighing more fast and stable.

ANALYTICAL BALANCE MAGNET SENSOR



MODEL	TL-200E	TL-300E	TL-500E	TL-600E	TL-300I	TL-500I	TL-600I
Capacity (g)	200g	300g	500g	600g	300g	500g	600g
Sensor	Magnetic						
Resolution (g)	0.001g (1mg)						
Min weighing (g)	0.004g						
Stable Time	≤3 Sec						
Display	LCD Display (White Back light and Black Font)						
Pan Size	Φ 80mm / Φ 90mm						
Operation temp	5 - 25 ℃						
Repeat ability	± 0.002g						
Liner	± 0.003g						
Draft shield size	240mm x 190mm x 165mm						
Cal.weight	External Calibration (200g)				Internal Calibration (Auto)		
Option	RS232C & RS485, Printer interface						
Packing Size	Inside:415x285x430mm			Outside:475x315x485mm			
G.weight	10 Kg						

