

ANALYTICAL PRECISION BALANCE

0.0001g

TOPLAB INDIA
"Serving Science with Solution"

NEW

TOPLAB INDIA'S Analytical Weighing Balances are High accurate weighing results, fast and stable and userfriendly operations. The Balance has Large LCD display white backlit brings clearer and more comfortable visual sense. The Balances has Electromagnetic sensor technology with External Calibraiton facility, Stainless Steel Pan, Glass shield.

Used in the Laboratories filed of Pharmaceutical, Chemical, Agriculture, Universities and Institutes, Scientific Research and other units for the rapid determination of quality and quantity of objects.

The Balance has Counting percentage(%) Overload alarm etc function, Full range tare, 18 kinds of units for choose, RS232 interface, Printer interface optional.

TECHNICAL SPECIFICATION:

MODEL	TL-220C
Capacity (g)	220g
Calibration	External (Internal Cal. Optional)
Resolution (g)	0.0001
Min Weighing (g)	0.0004
Stable Time	≤ 3 S
Display	LCD (white back light with black front)
Pan Size	Φ 80mm
Operation temp.	5-25℃
Repeatability	± 0.0002g
Linear	± 0.0003g
Draft shield size	240mmx190mmx265mm
Cal. weight	200g
Option	RS232 & RS485, Printer
Packing size	Inside:415x285x430mm Outside:475x315x485mm
G. weight	10kg

Note: To order internal calibration with Balance, Request for Model TL-220C-i

ANALYTICAL BALANCE MAGNET SENSOR

