

ANALYTICAL WEIGHING BALANCE

0.001g (1mg)

TOPLAB INDIA
"Serving Science with Solution"



TOPLAB INDIA'S Analytical Weighing Balances are High accurate weighing results, fast and stable and userfriendly operations. The Balance has Large LCD display white backlit brings clearer and more comfortable visual sense. The Balances has Electromagnetic sensor technology with External Calibraiton facility, Stainless Steel Pan, Glass shield.

Used in the Laboratories filed of Pharmaceutical, Chemical, Agriculture, Universities and Institutes, Scientific Research and other units for the rapid determination of quality and quantity of objects.

The Balance has Counting percentage(%) Overload alarm etc function, Full range tare, 18 kinds of units for choose, RS232 interface optional.

ANALYTICAL BALANCE



TECHNICAL SPECIFICATION :

Model No.	TL-200A	TL-300A	TL-500A
Max. Capacity(g)	200g	300g	500g
Technology	Electromagnetic Sensor		
Resolution(g)	0.001g (1mg)		
Min. Weighing(g)	0.010g		
Stable Time(Sec)	≤8s		
Display	Large LCD display white backlit		
Pan Size(mm)	φ80mm (Stainless Steel Pan)		
Operation Temp.	5°C - 25°C		
Repeatability(g)	±0.002g		
Linearity(g)	±0.003g		
Dimension(mm)	320x210x315 mm		
Cal. Weight(g)	200g		500g
Calibration	External Calibration		
Optional	RS232 interface		
G. Weight	3.6kg		
Power Supply	AC220V/50Hz		

