

Medical & Laboratory cold storage system

Medical & Pharmacy Refrigerator

Temperature:
2-8°C

Capacity:
50~1500L



Microprocessor control

Ensures the most accurate inside temperature at 2-8 °C and makes temperature resolution 0.1 °C.



Forced-air cooling

Reduces the physically induced temperate drop and enables an constant temperature to be maintained.



Door lock

Safety Lock to prevent unauthorized access which protects the storage in safe environment.




Heavy duty castors

Castors with break, easy moving and stop.



TOPLAB Accurate microprocessor control systems, excellent air cooling system and professional refrigerant system to ensure temperature 2~8°C Featuring Industrial design ensures heavy duty use and beautiful appearance. It is necessary equipment for hospital and laboratory to store drugs, vaccines, reagent, enzyme, etc.



TL-50V



Model	TL-50V	Capacity(L)	50 L
External Size(mm)	450 X 450 X 575	N.W./G.W.(kg)	14.5/15
Inner Size(mm)	382 X 355 X 415	Loading Quantity(20")	132
Package Size(mm)	530 X 520 X 630	Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Color sprayed steel
Power(W)	60	Inner Material	HIPS
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤46
Alarm System	High/low temperature, sensor failure		
Optional Accessories	Remote alarm system		

TL-60V



Model	TL-60V	Capacity(L)	60 L
External Size(mm)	470 X 450 X 670	N.W./G.W.(kg)	21/22
Inner Size(mm)	382 X 347 X 528	Loading Quantity(20")	96
Package Size(mm)	535 X 520 X 750	Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Color sprayed steel
Power(W)	60	Inner Material	HIPS
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤46
Alarm System	High/low temperature, sensor failure		
Optional Accessories	Remote alarm system		

TL-130V



Model	TL-130V	Capacity(L)	130 L
External Size(mm)	500 X 503 X 1292	N.W./G.W.(kg)	53/57
Inner Size(mm)	418 X 390 X 668	Loading Quantity(20")	36
Package Size(mm)	600 X 610 X 1400	Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Color sprayed steel
Power(W)	150	Inner Material	Embossed Al board
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤49.8
Alarm System	High/low temperature, sensor failure, door ajar, power failure(8h)		
Optional Accessories	Remote alarm system, temperature printer		

TL-150V



Model	TL-150V	Capacity(L)	150 L
External Size(mm)	595 X 585 X 837	N.W./G.W.(kg)	50/55
Inner Size(mm)	525 X 426 X 730	Loading Quantity(20")	
Package Size(mm)		Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Color sprayed steel
Power(W)	65	Inner Material	HIPS
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤48
Alarm System	High/low temperature, sensor failure, door ajar, power failure(8h)		
Optional Accessories	Remote alarm system		

TL-236V



Model	TL-236V	Capacity(L)	236 L
External Size(mm)	580 X 560 X 1785	N.W./G.W.(kg)	63/65
Inner Size(mm)	500 X 420 X 1090	Loading Quantity(20")	24
Package Size(mm)	710 X 665 X 1900	Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Color sprayed steel
Power(W)	254	Inner Material	Color sprayed steel
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤46
Alarm System	High/low temperature, sensor failure, door ajar, power failure(8h)		
Optional Accessories	Remote alarm system, temperature printer		

TL-300V



Model	TL-300V	Capacity(L)	300 L
External Size(mm)	620 X 560 X 1970	N.W./G.W.(kg)	72/75
Inner Size(mm)	580 X 420 X 1270	Loading Quantity(20")	24
Package Size(mm)	735 X 685 X 2075	Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Color sprayed steel
Power(W)	260	Inner Material	Color sprayed steel
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤46
Alarm System	High/low temperature, sensor failure		
Optional Accessories	Remote alarm system, temperature printer		

TL-656V



Model	TL-656V	Capacity(L)	656 L
External Size(mm)	1220 X 642 X 1885	N.W./G.W.(kg)	158/188
Inner Size(mm)	1100 X 454 X 1325	Loading Quantity(20")	13
Package Size(mm)	1290 X 695 X 1820	Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Color sprayed steel
Power(W)	443	Inner Material	Color sprayed steel
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤52
Alarm System	High/low temperature, sensor failure		
Optional Accessories	Remote alarm system, temperature printer		

TL-1006V



Model	TL-1006V	Capacity(L)	1006
External Size(mm)	1220 X 872 X 1885	N.W./G.W.(kg)	207/237
Inner Size(mm)	1100 X 684 X 1325	Loading Quantity(20")	10
Package Size(mm)	1290 X 925 X 1820	Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Color sprayed steel
Power(W)	472	Inner Material	Color sprayed steel
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤53
Alarm System	High/low temperature, sensor failure, door ajar, power failure(8h)		
Optional Accessories	Remote alarm system, temperature printer		

TL-1500V



Model	TL-1500V	Capacity(L)	1500 L
External Size(mm)	1800 X 750 X 1965	N.W./G.W.(kg)	190/217
Inner Size(mm)	1680 X 635 X 1412	Loading Quantity(20")	6
Package Size(mm)	1940 X 880 X 2035	Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Stainless steel
Power(W)	650	Inner Material	Stainless steel
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤53
Alarm System	High/low temperature, sensor failure, door ajar, power failure(8h)		
Optional Accessories	Remote alarm system, temperature printer		

TL-100C



Model	TL-100C	Capacity(L)	100 L
External Size(mm)	570 X 580 X 880	N.W./G.W.(kg)	30/32
Inner Size(mm)	450 X 375 X 705	Loading Quantity(20")	54
Package Size(mm)	665 X 615 X 920	Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Color sprayed steel
Power(W)	71	Inner Material	Embossed Al board
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤46
Alarm System	High/low temperature, sensor failure, power failure(48h)		
Optional Accessories	Remote alarm system		

TL-200C



Model	TL-200C	Capacity(L)	200 L
External Size(mm)	865 X 660 X 875	N.W./G.W.(kg)	35/38
Inner Size(mm)	735 X 440 X 685	Loading Quantity(20")	36
Package Size(mm)	955 X 700 X 935	Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Color sprayed steel
Power(W)	89	Inner Material	Embossed Al board
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤46
Alarm System	High/low temperature, sensor failure, power failure(48h)		
Optional Accessories	Remote alarm system		

TL-300C



Model	TL-300C	Capacity(L)	656 L
External Size(mm)	1120 X 745 X 875	N.W./G.W.(kg)	49/52
Inner Size(mm)	972 X 490 X 660	Loading Quantity(20")	22
Package Size(mm)	1205 X 790 X 935	Insulation	PURF
Voltage	220/110V; 50/60Hz	External Material	Color sprayed steel
Power(W)	110	Inner Material	Embossed Al board
Setting Temp.(°C)	2-8	Compressor/QTY	Brand/1
Ambient Temp. (°C)	10-32	Noise(db)	≤46
Alarm System	High/low temperature, sensor failure, power failure(48h)		
Optional Accessories	Remote alarm system		

Product by:

TOPLAB INDIA

Address:

601, 6th FLOOR, SHAH PRIMA,
PLOT NO: 13, SECTOR-02, KHARGHAR,
NAVI MUMBAI- 410 210, M.S., INDIA

E-mail: info@toplabindia.com

Website: www.toplabindia.com

Distributor: