

High Capacity Weighing Scales

(10kg/ 15kg / 20kg / 30kg / 50kg / 60kg / 100kg / 150kg / 300kg)

TOPLAB INDIA
"Serving Science with Solution"

WEIGHING SCALE

TOPLAB INDIA'S TL-xxxJS/TS Series Weighing Balances are manufactured by using high precision weighing sensors to keep weighing stable. It is very accurate and precision Balance in Industry. It is used for high capacity weighing in Laboratory, Factory, Ware house, Stores and many more industries.

It supports different interface and accessories to extend functions and the under weighing function is suitable for density tests.

Model	TL-10JS	TL-15JS	TL-20JS	TL-30JS	TL-50JS	TL-60TS	TL-100TS	TL-100TS	TL-300TS
Max. Capacity	10kg	15kg	20kg	30kg	50kg	60kg	100kg	150kg	300kg
Sensor Type	Load Cell								
Resolution (g)	0.1g	0.1g	0.1g	0.1g	1.0g	2.0g	5g	10g	20g
Min. Weighing (g)	5.0g	5.0g	10.0g	10.0g	200g	400g	400g	1kg	2kg
Stable Time:	≤ 2 Sec								
Display	Digital LCD Display								
Pan Size	280 x 210 mm				320 x 420 mm			400 x 500mm	
Operation Temp.	5 - 25℃								
Repeatability	± 2.0g				± 20g				
Linearity	± 3.0g				± 30g				
Column Height	N.A.				500mm			600mm	
Interface	RS232 interface (Optional)								
Dimensions	320x290x115mm				N.A.				
Weight	3.0kg				10.0kg				

*The calibration weights are optional available at extra cost on customer demand.

