

# ELECTRONIC WEIGHING BALANCE

## 0.01g

**TOPLAB INDIA**  
"Serving Science with Solution"

**NEW**

TOPLAB INDIA'S Analytical Weighing Balances are High accurate weighing results, fast and stable and userfriendly operations. The Balance has Large LCD display white backlit brings clearer and more comfortable visual sense. The Balances has Electromagnetic sensor technology with External Calibraiton facility, Stainless Steel Pan, Glass shield.

Used in the Laboratories filed of Pharmaceutical, Chemical, Agriculture, Universities and Institutes, Scientific Research and other units for the rapid determination of quality and quantity of objects.

The Balance has Counting percentage(%) Overload alarm etc function, Full range tare, 18 kinds of units for choose, RS232 interface, Printer interface optional.

### TECHNICAL SPECIFICATION :

Model No.	TL-500D	TL-1000D	TL-2000D	TL-3000D	TL-5000D
Max. Capacity(g)	500g	1000g	2000g	3000g	5000g
Technology	Electromagnetic Sensor				
Resolution(g)	0.01g				
Min. Weighing(g)	0.040g				
Stable Time(Sec)	≤8s				
Display	Large LCD display white backlit				
Pan Size(mm)	φ130mm				173X188 mm
Operation Temp.	5°C - 25°C				
Repeatability(g)	±0.02g				
Linearity(g)	±0.03g				
Dimension(mm)	320x210x175				320x210x95
Cal. Weight(g)	500g		1000g		2000g
Calibration	External Calibration				
Optional	RS232 interface ( optional)				
G. Weight	4 Kg				

### ELECTRONIC BALANCE

